

Clinical Performance of Reagent by Visual Analysis

CONTENTS

Reagent Pad Name	Page
Albumin	2
Bilirubin	3
Blood	4
Creatinine	5
Glucose	6
Ketone	7
Leukocyte	8
Nitrite	9
pH	10
Protein	11 - 12
Specific Gravity	13
Urobilinogen	14 - 15

Clinical Performance of Reagent by Visual Analysis

Reagent Pad	<i>Albumin (Protein-Low)</i>
Technical Name	<i>Albumin (1)</i>
Clinical Study Report	<i>CTD 99-06 *</i>
Date of Report	<i>June 1999</i>
Type of Study	<i>External, 2 clinical sites</i>
Test Solutions	<i>Clinical urines</i>
Number of Reagent Lots Used	<i>2</i>
Number of Results Obtained	<i>734</i>
Comparative Method	<i>Nephelometric immunoassay</i>
Maximum Concentration considered "Negative"	<i>< 20 mg/dL albumin</i>
Analytical Sensitivity	<i>69 %</i>
Analytical Specificity	<i>92 %</i>

* Located in Central File, Clinical Trials Department, Elkhart, Indiana

Clinical Performance of Reagent by Visual Analysis

Reagent Pad	<i>Bilirubin</i>
Technical Name	<i>Bilirubin (2)</i>
Clinical Study Report	<i>MPFO80685 *</i>
Date of Report	<i>August 1985</i>
Type of Study	<i>In-house evaluation</i>
Test Solutions	<i>Contrived urines</i>
Number of Reagent Lots Used	<i>2</i>
Number of Results Obtained	<i>144</i>
Comparative Method	<i>Urines at 4 fixed levels using assayed bilirubin stock standard</i>
Maximum Concentration considered "Negative"	<i>< 0.20 mg/dL bilirubin</i>
Analytical Sensitivity	<i>99 %</i>
Analytical Specificity	<i>89 %</i>

* Located in Central File, Clinical Trials Department, Elkhart, Indiana

Clinical Performance of Reagent by Visual Analysis

Reagent Pad	<i>Blood</i>
Technical Name	<i>Occult Blood (9)</i>
Clinical Study Report	<i>Technical Review, QAD 90-36 *</i>
Date of Report	<i>1990</i>
Type of Study	<i>In-house evaluation</i>
Test Solutions	<i>Contrived urines</i>
Number of Reagent Lots Used	<i>5</i>
Number of Results Obtained	<i>280</i>
Comparative Method	<i>Urine at 7 fixed levels using assayed blood stock standard</i>
Percent Agreement within \pm one color block	<i>100 % at all 7 Levels</i>

* Located in Central File, Clinical Trials Department, Elkhart, Indiana

Clinical Performance of Reagent by Visual Analysis

Reagent Pad	<i>Creatinine</i>	
Technical Name	<i>Creatinine (2)</i>	
Clinical Study Report	<i>CTD 99-06 *</i>	
Date of Report	<i>June 1999</i>	
Type of Study	<i>External, 4 clinical sites</i>	
Test Solutions	<i>Clinical urines</i>	
Number of Reagent Lots Used	<i>2</i>	
Number of Results Obtained	<i>1585</i>	
Comparative Method	<i>Jaffe Reaction</i>	
Percent Agreement within \pm one color block	<u><i>Color Block</i></u>	<u><i>% Agreement</i></u>
	<i>10 mg/dL</i>	<i>96</i>
	<i>50 mg/dL</i>	<i>97</i>
	<i>100 mg/dL</i>	<i>96</i>
	<i>200 mg/dL</i>	<i>99</i>
	<i>300 mg/dL</i>	<i>99</i>

* Located in Central File, Clinical Trials Department, Elkhart, Indiana

Clinical Performance of Reagent by Visual Analysis

Reagent Pad	<i>Glucose</i>
Technical Name	<i>Quantitative Glucose (1)</i>
Clinical Study Report	<i>MPFO80685 *</i>
Date of Report	<i>August 1985</i>
Type of Study	<i>External, 2 clinical sites</i>
Test Solutions	<i>Clinical urines</i>
Number of Reagent Lots Used	<i>2</i>
Number of Results Obtained	<i>2577</i>
Comparative Method	<i>Hexokinase Assay</i>
Maximum Concentration considered "Negative"	<i>49.9 mg/dL glucose</i>
Analytical Sensitivity	<i>99%</i>
Analytical Specificity	<i>98%</i>

* Located in Central File, Clinical Trials Department, Elkhart, Indiana

Clinical Performance of Reagent by Visual Analysis

Reagent Pad	<i>Ketone</i>
Technical Name	<i>Ketone (3)</i>
Clinical Study Report	<i>MPFO80685 *</i>
Date of Report	<i>August 1985</i>
Type of Study	<i>External, 2 clinical sites</i>
Test Solutions	<i>Clinical urines</i>
Number of Reagent Lots Used	<i>2</i>
Number of Results Obtained	<i>2576</i>
Comparative Method	<i>Sodium Nitroprusside Assay</i>
Maximum Concentration considered "Negative"	<i>9.9 mg/dL ketone</i>
Analytical Sensitivity	<i>92 %</i>
Analytical Specificity	<i>99 %</i>

* Located in Central File, Clinical Trials Department, Elkhart, Indiana

Clinical Performance of Reagent by Visual Analysis

Reagent Pad	<i>Leukocyte</i>
Technical Name	<i>WBC (3)</i>
Clinical Study Report	<i>MPFO80685 *</i>
Date of Report	<i>August 1985</i>
Type of Study	<i>External, 3 clinical sites</i>
Test Solutions	<i>Clinical urines</i>
Number of Reagent Lots Used	<i>2</i>
Number of Results Obtained	<i>4824</i>
<u>Comparative Method #1</u>	<i>Chamber Count</i>
Maximum Concentration considered "Negative"	<i>< 10 whole blood cells / μL</i>
Analytical Sensitivity	<i>90 %</i>
Analytical Specificity	<i>71 %</i>
<u>Comparative Method #2</u>	<i>Sediment Count</i>
Maximum Concentration considered "Negative"	<i>< 5 whole blood cells / HPF</i>
Analytical Sensitivity	<i>91 %</i>
Analytical Specificity	<i>69 %</i>

* Located in Central File, Clinical Trials Department, Elkhart, Indiana

Clinical Performance of Reagent by Visual Analysis

Reagent Pad	<i>Nitrite</i>
Technical Name	<i>Nitrite (4)</i>
Clinical Study Report	<i>MPFO80685 *</i>
Date of Report	<i>August 1985</i>
Type of Study	<i>External, 1 clinical site</i>
Test Solutions	<i>Clinical urines</i>
Number of Reagent Lots Used	<i>2</i>
Number of Results Obtained	<i>304</i>
Comparative Method	<i>Sulfanilamide-HTBQ assay</i>
Maximum Concentration considered "Negative"	<i>0.03 mg/dL nitrite</i>
Analytical Sensitivity	<i>88 %</i>
Analytical Specificity	<i>90 %</i>

* Located in Central File, Clinical Trials Department, Elkhart, Indiana

Clinical Performance of Reagent by Visual Analysis

Reagent Pad	<i>pH</i>	
Technical Name	<i>pH (1)</i>	
Clinical Study Report	<i>MPFO80685 *</i>	
Date of Report	<i>August 1985</i>	
Type of Study	<i>External, 2 clinical sites</i>	
Test Solutions	<i>Clinical urines</i>	
Number of Reagent Lots Used	<i>2</i>	
Number of Results Obtained	<i>2584</i>	
Comparative Method	<i>pH Meter</i>	
	<u><i>Color Block</i></u>	<u><i>% Agreement</i></u>
	<i>5.0 pH</i>	<i>99</i>
	<i>6.0 pH</i>	<i>99</i>
Percent Agreement within \pm one color block	<i>6.5 pH</i>	<i>90</i>
	<i>7.0 pH</i>	<i>89</i>
	<i>7.5 pH</i>	<i>98</i>
	<i>8.0 pH</i>	<i>89</i>
	<i>8.5 pH</i>	<i>100</i>

* Located in Central File, Clinical Trials Department, Elkhart, Indiana

Clinical Performance of Reagent by Visual Analysis

Reagent	<i>Protein</i>
Technical Name and Suffix	<i>Protein (4)</i>
Clinical Study Report	<i>CTD 94-04 *</i>
Date of Report	<i>May 1994</i>
Type of Study	<i>External, 4 clinical sites</i>
Test Solutions	<i>Clinical urines</i>
Number of Reagent Lots Used	<i>2</i>
Number of Results Obtained	<i>390</i>
<u>Comparative Method #1</u>	<i>human serum albumin assay</i>
Maximum Concentration considered "Negative"	<i>< 4.0 mg/dL albumin</i>
Analytical Sensitivity	<i>99.5 %</i>
Analytical Specificity	<i>67.4 %</i>
<u>Comparative Method #2</u>	<i>total protein assay</i>
Maximum Concentration considered "Negative"	<i>< 15 mg/dL total protein</i>
Analytical Sensitivity	<i>99.4 %</i>
Analytical Specificity	<i>72.4 %</i>

* Located in Central File, Clinical Trials Department, Elkhart, Indiana

Clinical Performance of Reagent by Visual Analysis

Reagent	<i>Protein (Protein-High)</i>
Technical Name and Suffix	<i>Protein (4)</i>
Clinical Study Report	<i>CTD 99-06 *</i>
Date of Report	<i>June 1999</i>
Type of Study	<i>External, 2 clinical sites</i>
Test Solutions	<i>Clinical urines</i>
Number of Reagent Lots Used	<i>2</i>
Number of Results Obtained	<i>1711</i>
<u>Comparative Method #1</u>	<i>total protein (TP) assay and human serum albumin (HSA) assay</i>
Maximum Concentration considered "Negative"	<i>Dual criteria: < 25 mg/dL TP and < 10 mg/dL HSA</i>
Analytical Sensitivity	<i>93 %</i>
Analytical Specificity	<i>91 %</i>
<u>Comparative Method #2</u>	<i>total protein assay</i>
Maximum Concentration considered "Negative"	<i><25 mg/dL total protein</i>
Analytical Sensitivity	<i>92 %</i>
Analytical Specificity	<i>85 %</i>

* Located in Central File, Clinical Trials Department, Elkhart, Indiana

Clinical Performance of Reagent by Visual Analysis

Reagent Pad	<i>Specific Gravity</i>	
Technical Name	<i>Specific Gravity (3)</i>	
Clinical Study Report	<i>MPFO80685 *</i>	
Date of Report	<i>August 1985</i>	
Type of Study	<i>External, 2 clinical sites</i>	
Test Solutions	<i>Clinical urines</i>	
Number of Reagent Lots Used	<i>2</i>	
Number of Results Obtained	<i>2558</i>	
Comparative Method	<i>Total Solids Meter</i>	
	<u><i>Color Block</i></u>	<u><i>% Agreement</i></u>
	<i>1.005</i>	<i>98</i>
	<i>1.010</i>	<i>77</i>
Percent Agreement	<i>1.015</i>	<i>82</i>
within \pm one color block	<i>1.020</i>	<i>71</i>
	<i>1.025</i>	<i>79</i>
	<i>1.030</i>	<i>59</i>

* Located in Central File, Clinical Trials Department, Elkhart, Indiana

Clinical Performance of Reagent by Visual Analysis

Reagent Pad	<i>Urobilinogen</i>
Technical Name	<i>Urobilinogen (4)</i>
Clinical Study Report	<i>Technical Review S85-40 *</i>
Date of Report	<i>1985</i>
Type of Study	<i>In-house evaluation</i>
Test Solutions	<i>Clinical urines</i>
Number of Reagent Lots Used	<i>2</i>
Number of Results Obtained	<i>216</i>
<u>Comparative Method #1</u>	<i>Watson's Aldehyde Assay</i>
Maximum Concentration considered "Negative"	<i>< 2.0 mg/dL urobilinogen</i>
Analytical Sensitivity	<i>47 %</i>
Analytical Specificity	<i>98 %</i>
<u>Comparative Method #2</u>	<i>Uro3 test pad</i>
Maximum Concentration considered "Negative"	<i>< 2.0 mg/dL urobilinogen</i>
Analytical Sensitivity	<i>100 %</i>
Analytical Specificity	<i>100 %</i>

* Located in Central File, Clinical Trials Department, Elkhart, Indiana

Clinical Performance of Reagent by Visual Analysis

Reagent Pad	<i>Urobilinogen</i>
Technical Name	<i>Urobilinogen (4)</i>
Clinical Study Report	<i>MPFO80685 *</i>
Date of Report	<i>August 1985</i>
<u>Type of Study (#1)</u>	<i>External, 1 clinical site</i>
Test Solutions	<i>Clinical Urines</i>
Number of Reagent Lots Used	<i>2</i>
Number of Results Obtained	<i>2119</i>
Comparative Method	<i>Uro3 test pad</i>
Maximum Concentration considered "Negative"	<i>< 2.0 mg/dL urobilinogen</i>
Analytical Sensitivity	<i>84 %</i>
Analytical Specificity	<i>100 %</i>
<u>Type of Study (#2)</u>	<i>In-house Evaluation</i>
Test Solutions	<i>Contrived urines at 4 fixed levels using assayed urobilinogen stock</i>
Number of Reagent Lots Used	<i>3</i>
Number of Results Obtained	<i>144</i>
Comparative Method	<i>Uro3 test pad</i>
Maximum Concentration considered "Negative"	<i>< 2.0 mg/dL urobilinogen</i>
Percent Agreement	<i>Normal: 100% Abnormal: 100%</i>

* Located in Central File, Clinical Trials Department, Elkhart, Indiana